


LCD-1000/1000P

Rexxam



LCD Chart

Smart, Slim & Easy to use

LCD-1000/1000P




LCD Chart


LCD-1000/1000P

FEATURES

- High resolution 24 inch display panel with circular polarization(LCD-1000P)/ 23 inch display panel(LCD-1000).
- Various types of standard charts, *polarized charts for Binocular test, ETDRS chart, Contrast test, *3D test, Color Vision test, Maddox test(LED), etc.
- Synchronize with Digital Refractor DR-900.
- Random function , Program mode, Masks, etc.
- 4 kinds of channels for the remote control are available.



Synchronize with Digital Refractor DR-900



LCD-1000/1000P



DR-900

SPECIFICATIONS

Monitor	Type	*24 inch / 23 inch TFT LED Backlight
	Resolution	1920(H) x 1080(V) Pixels
Functions	Chart	1. Landolt 2. Tumbling E 3. Numbers 4. Letters 5. Child 6. Special
	Masking	Single / Horizontal Line / Vertical Line
	Filter	Red and Green / *Circular Polarized Filter
	Special function	Mirror mode / Contrast / White and Black reverse Random / Crowd
	Working Distance	1.0m~8.0m (Step 0.25m) **
	Control device	Wireless IR remote controller / Digital refractor DR-900
	Screen saver function	Screen off / Chart show (OFF, 5sec, 10sec, 30sec)
	Screen saver waiting time	5min / 10min / 30min
Specifications	Body	573(W) x 363(V) x 46(D) mm / 4.3 kg
	Remote controller	64(W) x 195(V) x 20(D) mm
	Power supply (Main unit)	DC 12V 3.4A
	Power	100~240V
	Power consumption	70VA

LCD-1000/1000P Standard Accessories

- Operation manual
- Power cord
- AC adaptor
- Battery
- Wall mount
- R & G glass
- *Pola glass
- Remote controller

*only LCD-1000P

**Some charts will appear incomplete due to difference between screen resolution and image size.

Design and specifications are subject to change without notice.

Manufacturer _____

Distributed by _____

Rexxam

Rexxam Co., Ltd.

958 Ikeuchi, Konan, Takamatsu-shi,
Kagawa 761-1494 Japan

Contact _____

Eye-care Instruments Sales Dept. Tokyo Office

2-4-2 Kandatsukasa-machi, Chiyoda-ku, Tokyo, 101-0048, Japan

TEL:81-3-6262-9471 FAX:81-3-6262-9472

E-mail: eye@rexxam.co.jp

URL: http://www.rexxam.co.jp